

UltifreshTM Reusable Mask



















Breathable

Anti-Bacterial Washable

Cost Saving

BFE≥99%

Eco-Friendly

High Color Fastness Dust Free Production

tion Safe To Wear

PRODUCT SUMMARY



	FIT		FITS FACE CONTOUR DESIGN
	HYDROPHOBIC/ WATER REPELLANT	\checkmark	YES
Ø	FABRIC TYPE	\checkmark	KNITTED FABRIC THAT GOES WITH FACE CONTOUR
	BREATHABLE & COMFORT	\checkmark	HIGH
-8-8-	EASY TO WASH	\checkmark	DRIES IN 6 - 10 HOURS
	COLOR FASTNESS	\checkmark	ABLE TO WITHSTAND 40 - 60 DEGREE WATER
	ADJUSTABLE LOOPS	\checkmark	YES
	BACTERIAL FILTRATION		99.9% HIGH BFE PROTECTION
Ŵ	SAFE TO WEAR	\checkmark	CERTIFIED BY SGS & INTERTEK
	ANTI - BACTERIAL	\checkmark	99.94% SUSTAINABLE AFTER 50 WASHES
	NATURAL FABRIC TOUCHING OUR SKIN	\checkmark	100% NATURAL MADE COTTON
	DUST FREE ENCLOSED PRODUCTION	\checkmark	YES
Ŷ	CUSTOMIZE OPTION	\checkmark	ABLE TO ADD LOGO / CHANGE STRING / CUSTOM PACKING

• Drying speed based on air dry in room temperature


Number of usage based on normal washing

• Outer mask fabric is hydrophobic & water repellant based on World Health Organization(WHO) mask guidelines

FACE MASK MARKER

Kids Face maskAdult Face maskImage: state state

FACE MASK GUIDELINES





Adjustable Ear Loop

Suitable for all sizes








Extreme Comfort




99.9% high BFE protection



ACTUAL PRODUCT



3X SUPERIOR LAYERS PROTECTION





MASK FILTRATION EFFICIENCY



OUTERMOST LAYER

Ultifresh technology fabric.

bacteria.

capable of eliminating 99.94%


CENTER LAYER

BFE 99% melt-blown fabric to filter out contaminants.

O INNERMOST LAYER

100% soft style cotton. Giving all users maximum comfort.

ULTIFRESH[™] REUSABLE MASKS FILTER 99 PERCENT OF PARTICLES WITH A MASS MEDIAN DIAMETER OF 0.3 μm.



WASHABLE & REUSABLE





Hand washable



Do not tumble dry











HOW TO USE





SPREAD OPEN THE MASK




2

INSERT THE FILTRATION SHEET FROM THE SIDE



TUCK THE EARLOOPS OVER YOUR EARS



ADJUST THE MASK TO FIT YOUR FACE COMFORTABLY

4

COMPARING

ULTIFRESH[™] REUSABLE MASK TO WHO RECOMENDED FACE MASK

ULTIFRESH[™] REUSABLE MASK

OUTERMOST LAYER

Function:

Ultifresh treatment

hydrophobic (repels water)

Limits external contamination/ droplets from penetrating through to the wearer's nose and mouth. Material: Anti bacteria fabric with

CENTER LAYER

Function:

Prevents germs from entering our nose and mouth. Material: Virus filtration sheets

INNERMOST LAYER Function:

Absorb the heat and liquid droplets from our body, leaving our skin dry and comfortable.

Material:

100% softstyle cotton, hydrophilic (absorbs water)

DESIGNED AND MANUFACTURED ACCORDING TO WORLD HEALTH ORGANIZATION (WHO) RECOMMENDED FABRICS

INNER LAYER ____

Material properties: Hydrophilic (absorbs water) Material examples: Cotton or mixed cotton Function: Absorbs liquid droplets

MIDDLE LAYER _

Material properties: Made from synthetic non-woven material that prevents water absorption Material examples: Polypropylene or cotton layers Function: Increase the filtration or absorbs liquid droplets

OUTER LAYER

Material properties: Hydrophobic (repels water) Material examples: Polypropylene, polyester or a mix of both Function:

Limits external contamination from penetrating through to the wearer's nose and mouth

Source: Ministry of Health and WHO

STANDARD PROTECTIVE MASK

SGS

PROFESSIONAL TEST REPORTS

SGS and Intertek Certified

nterter

-

SGS

ntorto

Theorem of

SGS

SGS

-





BFE 99%

Official Recognition of The Medical Standard

OUTERMOST

LAYER



Melt-blown fabric to filter contaminants

99% Bacterial Filtration Efficiency

DUST FREE BACTERIA FREE CLOSURE ROOM




AVAILABLE DESIGNS

Adult Mask




PRODUCT DISPLAY



CUSTOMER REVIEWS

More than 1000+ positive review

Ratings & Reviews of 🔴 ** [READY 🍫 STOCK] (Ultifresh Reusable 🕄 Washable Mask) [1 Reusable Mask + 1 Virus Filter Sheet ...



Product Reviews

Filter: All star

3 weeks ago

Sort: Relevance

by Sharif-Sharmizi S. 💡 Verified Purchase

Item arrived very quickly, I am pleased with the material and how easy it is to fit the filter into the mask. Because it is also not made out of a hard material, it means it is easier to breathe in and comfortable to wear. Recommended.



Color Family:1 Navy Blue Color Reusable Mask + 1 Virus Filter Sheet



"Yes, happy with the quality and pretty comfy. It fits well too!" - Elson, Nicole & Christy in the Steel Industries (JB)



"Have my son to wear the mask as it offers maximum protection & fits his face contour with its construction." - Shermaine (SG) Researcher




CUSTOMIZE DESIGNS

Past Projects

Customize Label




ΥΙDEOJET



CUSTOMIZE DESIGNS

Past Projects

Full Dye Sub

















CUSTOMIZE DESIGNS

Past Projects

Custom Print/Embroidery



CUSTOM MADE Open For Printing of Own Logo / Design





POSITIONING OF OWN LOGO OR DESIGN AT MASK EDGE, LEFT CORNER



WHERE TO BUY



SINGAPORE

READY STOCK AVAILABLE AT THE FOLLOWING STORES

guardian

Lazada mothercare Resilis



MALAYSIA

READY STOCK AVAILABLE AT THE FOLLOWING PHARMACIES

- **Alex Pharmacy**
- **Ascen Plus Pharmacy**
- **Avanti Pharmacy**
- **Big Pharmacy**
- **Careplus Pharmacy**
- **Caring Pharmacy**
- Farmasi Howe Cheang Medical Supply
- Farmasi Murni Marketing
- Farmasi Sungai Pelek **Healthcare Enterprise**

- First Pharmacy
- Georgetown Pharmacy Langkawi
- Georgetown Medical Supplies
- Gold Coast Pharmacy
- Gold River Pharmacy
- Green Trees Pharmacy
- Happy Pharmacy
- Healthway Pharmacy
- HTM Pharmacy

- IJ Pharmacy
- **Mediland Pharmacy**
- Mega Kulim Pharmacy
- MY Pharmacy
- Pharmacy Alley
- Pharmhouse Farmasi
- Pure Pharmacy
- SK Pharmacy
- TMC Pharmacy
- Wah Choi Pharmacy
- Wellings Pharmacy



ORDER FORM



Singapore Head Office

Address : 23 Genting Road, #03-03, Chevalier House, Singapore 349481 Operation Hours : Monday - Friday | 09:00a.m. to 06:00p.m. (Local Time) Phone No. : +65 6282 2858 | Fax : +65 6282 2856

Malaysia Branch Office

Address : No.15 Jalan Desa Tropika 1, Taman Perindustrain Tropika, 81800 Ulu Tiram, Johor, Malaysia Operation Hours : Monday - Friday | 08:30a.m. to 05:30p.m. (Local Time) Phone No. : +607 213 2078 | +6011 3178 1707

赣州优迪飞服饰有限公司

Address:中国江西省赣州市于都县楂林工业园龙门路北侧, 怡信大道东侧(怡信路57号4楼无尘车间)